



Quality Assurance of the Hyper-Kamiokande Outer Detector Photomultiplier Tubes



Robert Kralik

Hyper-Kamiokande Outer Detector


Hyper-Kamiokande Outer Detector




Hyper-Kamiokande (HK) detector



3" Outer Detector (OD) PMTs



Lined with reflective tyvek sheets



Outer Detector overview



Partially contained neutrino interaction

- OD near corners helps differentiate from corner-clipping muons

Cosmic muons , 45Hz = 4million/day

- Background to all physics measurements
- Ideally veto all, in practice best you can do without changing OD volume is $\sim 1:10^5$ level
- Fast OD veto avoids recording and reconstructing all events

JPARC beam neutrinos

- Far side provides separation of FC and PC events

JPARC beam neutrinos


- Sensitive to δ_{cp} & Δm_{3i}^2
- Near side provides separation for rock interactions

Passive shielding from radioactivity in rock

- Main value to Low-E analyses, does not require PMTs

Upwards muons from ν_μ interactions in rock

- Sensitive to Mass Ordering & Δm_{3i}^2
- Bottom endcap OD separates these from (lower-E) PC events




Outer Detector PMTs

Two candidate 3" Photomultiplier Tubes (PMTs) for the OD – Hamamatsu R14374 and NNVT N2031

- Tender now ongoing – decision in October 2025
- New PMT versions already in Kamioka – shipped to UK later
- OD will have 3600 3" PMTs in total
- We mostly care about single Photo Electron (PE) sensitivity
 - Most OD PMTs will receive <5 PEs for majority of muons
- Main parameters important for the OD PMTs:
 - Low dark rate (<2 kHz for 25°C)
 - High Quantum Efficiency (QE) across Cherenkov spectrum (both direct and reflected)
 - Good Collection Efficiency (CE) across the PMT including edges where WLS plates connect

More information in **TN0046** - "3-inch Photomultiplier Tubes for the Hyper-Kamiokande Outer Detector"




PMT Quality Assurance

Quality Assurance



Each PMT undergoes Quality Assurance
as it arrives
~450/month
PMTs are re-packaged and temporarily
stored afterwards

Transport to HK for construction
Start December 2026



- 2 Dark Boxes with each 14+1 PMTs each per day
- One reference PMT in each box used for validation
- Up to 140 PMTs/week → up to 560 PMTs/month
- Additional measurements require alternative setup
 - Considering measuring PMT+WLS plate setup, angular dep. of collection efficiency, multiple wavelengths

Measurements for each PMT

- Gain variation with input voltage
- Relative Quantum Efficiency (single wavelength or more?)
- Single Photo Electron characteristics (charge resolution, peak-to-valley ratio)
- Dark rate
- Maybe time resolution?





Quality Assurance Space

9 m

Spaces secured for OD at University of Toyama by logistics group

- Candidate (1F) already available: 8m x 9m
- Only need about 3x5 m for the QA
- Could also be used for assembly and temporary storage of PMTs
- Not attached to the storage area – requires transportation



Loads of extra space

Quality Assurance Setup


2 Dark boxes with 14 (+1 reference) PMTs each /day

- 3D printed PMT stands designed to support the cable rolls on the side
- Hamamatsu rolls measured to be $\varnothing \approx 30\text{cm}$
- Splitter boards inside the box
- Two separate LEDs for each dark box
- LED input split into 8 output fibres


Thorlabs 1-to-7 fibre bundle
BF74LS01



- Each LED output illuminates 4 PMTs
 - Not uniform but should be good enough
 - Also considering output from side of PMTs




Quality Assurance Setup



Input voltage and output signal from CAEN VME box

- 5x HV power supply (6 ports each) CAEN V6521P
- 2x ADC (16 ports each) CAEN VX1730SB
- VME to USB bridge CAEN VX3718



ADC and HV tested at KCL during PMT measurements

Quality Assurance procedure – Monday to Friday



Morning

- | | |
|---|------------|
| 1) Finish the dark rate analysis from the previous batch and send results to database | 15 min |
| 2) Unload previous batch of PMTs from the dark boxes and re-box them | 0.5 – 1 hr |
| • Move to storage space right away or wait for assembly at the QA space? | |
| 3) Unbox and load new batch of PMTs into dark boxes | 0.5 – 1 hr |
| • Label each PMT and its cable with a bar code | |
| 4) Turn on HV to initial value and leave them to warm up | 1 – 4 hr |

Afternoon:

- | | |
|---|------------|
| 1) Turn on LED and measure gain variation vs input voltage | 0.5 – 1 hr |
| 2) Analyze the gain variation and determine operation voltage (V_{op}) for each PMT | 10 min |
| 3) Set each PMT to this voltage and measure single PE characteristics and relative QE | 0.5 hr |
| 4) Turn off the LED and measure dark rate over night / over weekend | 10 hr |


Quality Assurance – other measurements



Considering additional measurements to test collection efficiency, WLS plates, multiple wavelengths...

- Would likely need to dedicate another dark box (and space) for this
- Pre-calibration of OD PMTs likely only towards the end of installation
 - Too late to discover any issues

- [Idea] For QA we will build a device similar to that on previous page: geomagnetic-field-compensated dark box with a PMT mount, optical switch, and light source
 - → Include check of PMT/WLS relative yield for *all* delivered PMTs ideally
 - Maybe only need to do this for first 5% of tubes if no problems
 - Will incur more measurement time than nominal QA if done one-by-one [needs thought]



- Do this test in the same location as nominal QA
- Still TBD, but Toyama U. Is preferred location

Quality Assurance – status



Previous KCL postdoc (Stephane) designed the setup, procured most of the equipment

- HV power supply and Analog-Digital Converters (ADCs) already used for PMT measurements
- Two working prototypes of PMT stands 3D printed and currently tested
- Several LEDs+drivers, fibres, fibre bundles, optical feedthroughs (various light injection setups)

To-do list

- 1) Finalize and test the setup
 - 3D printing design for the PMT support (specifics depends on choice of PMT), light injection design, PMT organization inside the dark box
- 2) Develop DAQ software, user interface, finalize analysis scripts, implement interface to HK database
- 3) Finish the procurement of all the coaxial cables, feed-throughs, splitter boards, DAQ computer

Testing and development (re)starting now until summer 2025

Ship and start validation on-site between end of 2025 and beginning of 2026

Analysis

PMT waveform and charge distributions



LED pulses trigger PMT readout

Collect 500000 waveforms; for each

- 1) Find minimum
- 2) Calculate and subtract baseline
- 3) Add bins in pre-defined signal range
- 4) Convert to total collected charge

$$Q \text{ [pC]} = \text{Current [mA]} * \text{time [ns]}$$

$$I \text{ [mA]} = \text{Voltage [mV]} / \text{Impedance } [\Omega]$$

Time ... pre-defined time window


Impedance = 50 Ω (at ADC)

Voltage... size of signal


Create histogram of integrated charges

Fit to get the desired characteristics



Single PE



Multi PE



↓ Integrate



Single PE Charge Distribution Fits

Using ROOT fit based on the Full PMT response function as in [https://doi.org/10.1016/0168-9002\(94\)90183-X](https://doi.org/10.1016/0168-9002(94)90183-X)

$$= \text{Poisson} \otimes (\text{Signal} * \text{Gaussian} + \text{Background} * (\text{Gaussian} \otimes \text{Exp}))$$

- Previously used a sum of Gaussians and an Exponential

$$\text{Gain} = (Q_1 - Q_0) / e$$


- Set to 5×10^6 by adapting input voltage

$$\text{Peak-to-Valley ratio} \equiv P/V$$

- Both Peak and Valley taken from histogram
- Requirement: $P/V > 2$

$$\text{PE resolution} = \sigma_1 / Q_1$$

- Previously used wrong equation:
 $\sigma_1 / (\text{Amp}_{1\text{PE}} * \sigma_1 * \sqrt{2\pi})$
- Requirement: $< 50\%$



Dark Rate and Quantum Efficiency calculation



Dark Rate is calculated as $DR = N_{\text{OverThreshold}} / (N_{\text{waveforms}} * \text{wavesize [s]})$

- $N_{\text{OverThreshold}}$ is the number of times all waveforms crossed below the threshold
- Wavesize is just the length of each waveform in seconds
- Error on the Dark Rate is calculated as $\sqrt{N_{\text{OverThreshold}}} / (N_{\text{waveforms}} * \text{wavesize [s]})$

Quantum Efficiency (QE) is calculate using a reference PMT: $QE_{\text{Measured}} = NPE_{\text{Measured}} * (QE_{\text{Reference}} / NPE_{\text{Reference}})$

- NPE is the calculate number of photo electrons
- In real life we really measure the detection efficiency (Quantum Efficiency * Collection Efficiency)

The analysis scripts are located in <https://github.com/rkralik5/PMTFitting>

- **PMTFit.C** = single PE fits
 - Already works very well
 - Currently finishing comparisons with Federico from RAL (gain calculation)
 - Possible to include time resolution fits
 - Can be used to get gain dependence on voltage
- **MultiPEFit.C**
 - Simple Gaussian fit to the multi-PE charge distribution
 - Finishing comparisons with Federico to ensure stability and reliability
 - Federico believes this is more reliable for gain calculation
- **DarkRate.C**
 - Should already work well enough

Setup design

List of equipment available [here](#)

- Two dark boxes - one assembled
 - Still need to get someone to machine holes for feedthroughs
- Prototypes of PMT stands – different for Hamamatsu and NNVT
 - Hamamatsu stand is currently horizontal
- Several PMTs (4* Hamamatsu, 1 NNVT) [link to spreadsheet](#)
 - Neither of them final design
- Some LEDs and optical fibres, feedthroughs, fibre bundles...
 - Stephane's idea explained in previous slides
 - Other options for light input: separate fibre for each PMT, [diffuser ball](#), [laser ball](#)
 - Do we want multiple wavelengths? Are the LEDs enough?
 - Fibre holder probably 3D printed
- We have NNVT splitter boards – final design will probably use official HK OD ones



- 1) Learn how to use the 3D printer and collect all the PMT stand design prototypes
- 2) Design a vertical PMT stand for Hamamatsu
- 3) Design a fibre holder
- 4) Test Stephane's idea of light input
- 5) Test other fibre inputs
- 6) Adapt the design to be simple for QA shifters
- 7) Adapt the design so the individual PMT stands fit neatly together

Data Acquisition (DAQ)

Not developed yet – using CAEN software for the PMT measurements

- CAEN Scope – simple waveform measurement with triggers
- CAEN COMPASS – much more complicated software allowing for analysis (not necessary)
- CAEN WaveDump – should be a very simple software for ADC readout

The idea is that we will develop our own simpler case-specific software that will do the same thing

- CAENVMELib – list of functions to interact with the VME bridge
- CAENComm Library - Interface library for CAEN Data Acquisition Modules
- CAENDigitizer Library - Library of functions for CAEN Digitizers high level management

For the High Voltage control, we can either use the default software

- CAEN GECO

There are other possible libraries for the power control, we need to explore them and find out what is best

Starting to-do list



- 1) Read on the QA from the technote and this presentation
- 2) Look at CAEN WaveDump (for VX1730SB ADC) and try installing it on MAC
- 3) Look at the CAEN Lib documentation
- 4) Look at the documentation for the ADC
- 5) Access the VME bridge with simple commands
- 6) Get a waveform from a PMT
 - 1) With WaveDump
 - 2) With CAEN commands
- 7) Make a python (or other) script that will run the waveform reading for you, set every parameter, and save the waveforms

Useful links



- [PMT basics by Photonis](#)
- [3-inch Photomultiplier Tubes for the Hyper-Kamiokande Outer Detector - HK TN0046](#)
- [CAEN VME VX1730SB documentation](#)
- [CAEN HV V6521 documentation](#)

Thank you!